

VFI-1012-OP-0004 V1

## Unboxing Instructions for Visionflex Telehealth Cart

Cart Type	Ergotron SV40 with Drawer
AIO Type	Lenovo AIO Computer 24"
ProEX-Hub	No
PTZ Camera on Pole	Yes
Patient facing monitor	Yes
Document Version	1

Your Telehealth cart will arrive in a custom cardboard box, prewired and preconfigured. Depending on the selected configuration of the cart, some assembly will be required to mount the AIO Computer Screen onto the cart. The instructions below will guide you through this process, step by step.

**Tools required:** Scissors to cut packing straps and a knife to open the boxes.

---

### 1. Opening the Cart box.

Inspect the packaging and report any damage to Visionflex

Cut the strapping with scissors and open the top of the cart box with a sharp knife.



VFI-1012-OP-0004 V1

---

**2. Removing packaging and accessories from the box.**

Open the Box and remove internal packaging material being careful to separate any accessories which may have been packed inside the cart box.



**3. Tilting the Cart box upright**

Once the packaging contents have been removed and all the accessories have been carefully separated to one side, the box can be tilted upright so that the cart can be wheeled out of the box.



**4. Removing the Cart from the box**

Wheel the cart from the box. The wheels may be locked so please check and unlock if necessary.



VFI-1012-OP-0004 V1

### 5. Checking the Cart once removed from the box.

Once the cart has been removed from the box, inspect for damage and report any findings to the Visionflex support team.

You will notice that the cart has been prewired.



### 6. Mounting the monitor

Remove the monitor from the transport box and place face down on a smooth surface. Remove the 4 mounting bolts from the rear side of the monitor.

**Note:** This task can be more easily carried out by two people.



### 7. Washers for screws

Make sure that each thumb screw has the associated washer fitted.



## 8. Inserting the screws


Hold the monitor against the silver VESA mount on the cart and secure in place with the 4 thumb screws.

**Note:** This task can be more easily carried out by two people.



## 9. Connecting the monitor

Connect the preconfigured cables to the rear-side of the AIO Computer. These will include the power cable and the USB feed from the USB-Hub on the cart.



## 10. Connecting the Patient Facing Display

Install the patient facing display monitor that has a pre-installed mount on the back with the Swing Arm on the Cart.

Screw 2x Hexagonal screws on the top and 1x thumb screw on the bottom. Use provided hexagonal screw Wrench spanner. Plug power cable and HDMI cable on it.



### 11. Connecting PTZ camera.

Unscrew thumb screw (Green plastic circular level) of the ball head camera mount on the telescopic pole. Install PTZ camera on the top of the Telescopic pole Bracket. PTZ camera is mounted already with the base of Ball head camera Mount. Adjust the camera and screw tightly the thumb screw (Green plastic circular level). Plug power cable and USB cable on the back of the PTZ camera.



### 12. System Test


Connect the cart to mains power and turn on the AIO by pressing the power button on the right side of the device.

The username for the AIO computer is set to “**Telehealth Cart**” by default with no PW but your IT department can change this as required.



### 13. Congratulations

The support team at Visionflex look forward to working with you on the implementation of this new technology.



For more information or assistance with your cart assembly or to book in a training session, please contact the support team by email on [support@visionflex.com.au](mailto:support@visionflex.com.au) or by phone on 02 8914 4000.