

# EAR INFRARED THERMOMETER

## TD-1261C/D

### FEATURES

- Fast ear scan in one second
- Dual-color backlight
- Probe cover ejection and detection
- Memory capacity up to 10 results
- Bluetooth data transmission function (TD-1261C)
- Supplied with 16 ear covers!



### THERMOMETER SPECIFICATIONS

Temperature Measurement Sites	Ear Temperature
Temperature Range	32°C - 43°C
Accuracy	±0.2°C (±0.4°F) for the measurement at 36°C-39°C (96.8°F-102.2°F) ±0.3°C (±0.5°F) outside above temperature range
Display Resolution	0.1°C / 0.1°F
Temperature Unit	°C or °F set by factory
Backlight Screen	32°C 37.9°C 43°C Lo Normal Fever Hi
Probe Cover Detector	Equipped
Memory Capacity	10
Power Saving	3 min
Power Source	2 x AAA
Battery Life	1000 measurements
Dimension	164.8 (L) x 32.4 (W) x 38.7 (H) mm
Data Transmission	Bluetooth smart (TD-1261C); No Data Transmission (TD-1261D)
Weight	85.5g



www.taidoc.com

