


All-in-One Telehealth Cart


Your complete telehealth solution.

The AIO telehealth cart solution features a ProEX Desktop, a AIO medical grade computer, and an ergonomic, easy-to-manoeuvre medical cart.

Designed to meet multiple purposes, this universal cart solution combines the functionality of an all-in-one computer with the clinical capabilities of the ProEX Desktop.



ERGONOMIC



MULTI-PURPOSE



FLEXIBLE

Ergonomic, flexible,
and multi-purpose



All-in-One Telehealth Cart

Your complete telehealth solution.

The telehealth cart is an innovative, ergonomic and easy-to-manoeuvre medical cart that features sit-to-stand adjustment, independent screen positioning, side-to-side keyboard motion and an adjustable back-tilt keyboard tray.

The cart includes secure CPU compartment, scanner holder bracket, mouse holder, wrist rest, left or right mouse tray, rear handle, storage basket, cable storage and four lockable casters.

THE AIO COMPUTER

The medical grade AIO computer integrates a 24" screen with an energy efficient 10th gen Intel® Core™ i5 platform in a fan-less, slim, bacteria-resistant enclosure.

Designed for medical and healthcare applications, this LCD-integrated system provides space-saving solutions for health professionals to monitor, record and retrieve patient information as well as other point-of-care applications.

SPECIFICATIONS

Wheel-to-wheel base dimensions	W430 mm x D533 mm
Workspace dimensions	W567 mm x D502 mm
Shipping dimensions	W567 mm x D502 mm
Weight	26.8kg - unladen (39.1kg shipping weight)
Maximum load	15.9kg
Mounting standards	ESA FDMI MIS-D, 100/75, C (hole pattern: 100 x 100 mm & 75 x 75 mm)



For more information or to book a demonstration,
please contact Visionflex at sales@visionflex.com

+61 2 8914 4000

 visionflex.com



visionflex
Equality healthcare